

Sets Composition Configuration (5 - 2 - 2 - 1)

Number of weights									
12	20	21	22	24	25	26	28	29	30
Set from 1mg to									
500 mg	50 g	100 g	200 g	500 g	1 kg	2 kg	5 kg	10 kg	20 kg
1mg	1mg	1mg	1mg	1mg	1mg	1mg	1mg	1mg	1mg
2	2	2	2	2	2	2	2	2	2
2	2	2	2	2	2	2	2	2	2
5	5	5	5	5	5	5	5	5	5
10	10	10	10	10	10	10	10	10	10
20	20	20	20	20	20	20	20	20	20
20	20	20	20	20	20	20	20	20	20
50	50	50	50	50	50	50	50	50	50
100	100	100	100	100	100	100	100	100	100
200	200	200	200	200	200	200	200	200	200
200	200	200	200	200	200	200	200	200	200
500	500	500	500	500	500	500	500	500	500
	1g	1g	1g	1g	1g	1g	1g	1g	1g
	2	2	2	2	2	2	2	2	2
	2	2	2	2	2	2	2	2	2
	5	5	5	5	5	5	5	5	5
	10	10	10	10	10	10	10	10	10
	20	20	20	20	20	20	20	20	20
	20	20	20	20	20	20	20	20	20
	50	50	50	50	50	50	50	50	50
		100	100	100	100	100	100	100	100
			200	200	200	200	200	200	200
				200	200	200	200	200	200
					500	500	500	500	500
						1 kg	1 kg	1 kg	1 kg
							2	2	2
								2	2
							5	5	
								10	
								10	
								20	
Set from 1g to									
50 g	100 g	200 g	500 g	1 kg	2 kg	5 kg	10 kg	20 kg	
Number of weights									
8	9	10	12	13	14	16	17	18	

INSTRUMENTATION SERVICES, Inc.

281 Piñero Avenue, Río Piedras, PR 00927. Tel: (787) 765-5564

Fax: (787) 765-5585 e-mail: info@insco.us

INSCO METROLOGY, INC.

10434 NW 31st. Terrace, Miami, FL 33172 Tel. (305) 994-8031

Fax: (305) 994-7365 e-mail: miami@insco.us

INSCO DE MEXICO

Bldv. Toluca No. 43-C, CP53500 Naucalpana de Juarez

Edo. de México, México.

Tel: (525)55-359-0088 Fax: (525)55-358-3913 e-mail: info@inscomex.com